


Screen-Printed Co-Phthalocyanine/Carbon Electrode

Ref. 410



Disposable **Screen-Printed Co-Phthalocyanine/Carbon Electrodes (ref. 410)** are ideal for the determination of **hydrogen peroxide at a low detection potential**. These electrodes are recommended for the development of **enzymatic biosensors based on oxidases**, for working with microvolumes and for decentralized assays.

Ceramic substrate: L33 x W10 x H0.5 mm

Electric contacts: Silver

The electrochemical cell consists on:

Working electrode: Cobalt-Phthalocyanine/Carbon (4 mm diameter)


Auxiliary electrode: Carbon

Reference electrode: Silver


Co-Phthalocyanine/Carbon Electrodes are commercialised in a 75 units pack. They should be stored at room temperature, protected from light in a dry place.

Screen-Printed Co-Phthalocyanine/Carbon Electrode

Ref. 410



Amperometric detection of hydrogen peroxide in a flow injection analysis system with our easy to use flow cell (ref. FLWCL). The amperometric responses for decreasing H_2O_2 concentrations at a ref. 410 electrode show neither fouling nor memory effects. $E_{det} +0.4$ V; Flow rate 2.2 ml/min; Flow carrier 0.1 M phosphate buffer, pH 7.2



Calibration curve for hydrogen peroxide (in a 0.1 M phosphate buffer pH 7.2) from $1 \cdot 10^{-6}$ M to $2.5 \cdot 10^{-4}$ M in a FIA system (ref. FIASYSYSTEM using our Screen-printed Co-Phthalocyanine/Carbon electrodes. $E_{det} +0.4$ V; Flow rate 2.2 ml/min

Also, specific connectors that act as an interface between the screen-printed electrode and any potentiostat (refs. DSC, CAC) and other accessories are available at [DropSens](https://www.dropsens.com).

Related products



710



DSC



FLWCL



CELL



STAT400



STAT8000

Full Catalogue



Parque Tecnológico de Asturias - Edif. CEEI. 33428 LLanera (Asturias). Spain
(+34) 985 27 76 85 - info@dropsens.com - www.dropsens.com

Contact Form

