

# Eight channel boxed connector for $\mu$ Stat-i MultiX

## Ref. 4MMHCAST8

01

Eight channel boxed connector suitable to improve the set-up configuration acting as an interface between  $\mu$ Stat-i MultiX and screen-printed electrodes.

On one side, it has eight groups of 4 mm female banana connectors (WE, WE2, RE, AE and GND) to connect directly to the I-CABSTAT cable terminals.

On the other side it has eight five-pin connections recommended to link CASTDIR (SPE connector), BICASTDIR (dual SPE connector) or any other Metrohm DropSens cable terminated in a five-pin connector e.g. CAST-TLFL.



Dimensions: 30.0 x 10.0 x 5.0 cm (L x W x H)

[www.metrohm-dropsens.com](http://www.metrohm-dropsens.com)