

μ Stat 8000/P - μ Stat 4000/P Cable Connectors for any kind of electrodes

Refs. CABSTAT1
CABSTATMULTI

Flexible cables 1 m long which operate as an interface between μ STAT 8000/P or μ STAT 4000/P and any kind of conventional electrode, by means of 2 mm bananas and alligator clip connections.

μ Stat 8000/P - μ Stat 4000/P One Node Cable Connector for electrodes

Ref. CABSTAT1

Connects one node of μ STAT 8000/P or μ STAT 4000/P to any kind of electrochemical cell formed by one working, one auxiliary and one reference electrode.



μ Stat 8000/P Multichannel Cable Connector for electrodes

Ref. CABSTATMULTI

Connects μ STAT 8000/P Multichannel to any kind of electrochemical cell formed by up to eight working electrodes, one auxiliary and one reference electrode.



Also available for Ref. μ STAT 4000/P Multichannel under ref. CABSTATMULTI4.

Related products



STAT4000



STAT8000



4W110



STAT8000

Full Catalogue



Parque Tecnológico de Asturias - Edif. CEEI. 33428 LLanera (Asturias). Spain
(+34) 985 27 76 85 - info@dropsens.com - www.dropsens.com

Contact Form

