

Boxed connector for screen-printed electrodes Boxed connector for dual screen-printed electrodes

Ref. **DSC**, **BIDSC**, **DSC-P**, **DSC4MM**, **BIDSC-FET**



Dimensions 8.0 x 4.0 x 3.5 cm. (Length x Width x Height).

It operates as an interface between DropSens screen-printed electrodes (SPEs) and interdigitated electrodes with any kind of potentiostat.

Recommended when working by placing microvolume drops on our SPEs.

- Ref. **DSC** connects single (1 WE) SPEs
- Ref. **BIDSC** connects dual (2 WE) SPEs
- Ref. **DSC-P** connects single (1WE) plastic substrate SPEs
- Ref. **DSC4MM** connector with 4mm banana / crocodiles connects single (1WE) SPEs
- Ref. **BIDSC-FET** connects AUFET electrodes

www.metrohm-dropsens.com