



TEM image of Mesoporous Carbon Solution

DropSens Mesoporous Carbon Solution is homogeneous, black liquid with a high purity mesoporous carbon nanopowder. Carbon nanoparticles have a surface area $>200 \text{ m}^2/\text{g}$ and 64 \AA of average pore diameter.

The product is available in 1.0 or 5.0 mL order quantities. It has been specially designed for electrodes modification.

Related products



110MC



CNTSOL



GPHSOL



GPHXSOL

Full Catalogue



Parque Tecnológico de Asturias - Edif. CEEI. 33428 LLanera (Asturias). Spain
(+34) 985 27 76 85 - info@dropsens.com - www.dropsens.com

Contact Form

