


Membrane for Screen-Printed Electrodes


Ref. MEMB



This mesh is designed to perform **analysis** with **small volumes** of solution or sample (7,5-15 μ L), allowing a good covering of the 3 electrodes system in our Screen-Printed Electrodes.

Just place it over your SPE and work as usual but placing small volumes of sample.

The woven open mesh is a **monofilament fabric of polyamide**.



a_o : 35%

This membrane fits with the most of our screen-printed electrodes in catalogue. It's especially suitable for improving the wettability of the three electrode system, working with small volumes or delimiting the electrochemical cell.

Related products



110



220AT



BI10



P10



110CNT

Full Catalogue



Parque Tecnológico de Asturias - Edif. CEEI. 33428 LLanera (Asturias). Spain
(+34) 985 27 76 85 - info@dropsens.com - www.dropsens.com

Contact Form

