

Peristaltic Pump

Ref. PPUMP



DropSens Peristaltic Pump has a 12 roller pump head for 4 liquid flow channels adjustable from 0.005 ml/min to 30 ml/min. It is very easy to use, since pump tubing is placed around the pump head. Up to 4 tubes with similar or different diameters can be used at the same time, so up to 4 liquids can be pumped at similar or different flow rates.

An ON/OFF switch is available in the rear panel of the pump as well as the power connector to the means. In the front panel a rotational controller allows to select and modify easily the speed of the pump head (rpm), right or left rotation can be also selected in the front panel and a high speed button allows to work at maximum speed during a short period of time. Rollers configuration in the pump head allows an extremely reliable adjustment of liquid flow rate.

Sets of pump tubing with different diameters are available.

Dimensions: 160x160x200 mm

Weight: 4.5 kg

Related products



STAT400



STAT200



FLWCL



FIAEC

Full Catalogue



Parque Tecnológico de Asturias - Edif. CEEI. 33428 LLanera (Asturias). Spain
(+34) 985 27 76 85 - info@dropsens.com - www.dropsens.com

Contact Form

