

Product use and handling instructions

SPCEs STR

Please read these instructions carefully prior to beginning the assay:

The screen-printed carbon electrodes (SPCEs) modified with Streptavidin packaged in an aluminium envelope, heat sealed at the four corners, with a notch on each end of the top. To open the envelope and remove the electrodes properly, hold the envelope by its ends and break through the notches.

Before using the SPCEs modified with Streptavidin, the electrode surface should be washed with a buffer determined by the end user.

The electrodes should be stored at 2-8 °C.

Warning:

Do not hold the envelope through the center. It could damage the modification with Streptavidin.